Fixcet[®] Spinal Facet Screw System





- > Available in single-thread and dual-thread screws
- Cannulated screws for percutaneous approach
- Different thread pitches provide increased superior to inferior facet compression

Single-Thread Screw

Available in 25 - 45mm lengths (5mm Increments)

Dual-Thread Screw

Available in 30 - 45mm lengths (5mm Increments)

- 4.5mm diameter screws
- All screws cannulated
- 12mm percutaneous incision
- No washer required
- Simple instrumentation
- Various lengths for different spinal anatomies
- Lower profile than traditional pedicle screws
- Only 2 small incisions required for bi-lateral fixation
- Graft packer included for easy graft insertion

Different thread pitches provide increased superior to inferior facet compression

> Screw heads used for screwdriver engagement

Washer not required thus reduces profile

Anti-migration indent

Top threads provide additional locking without need of a washer

Cannulated screws for percutaneous approach

Optional 10mm washer is also available

Fixcet[®] Implants

Item Number	Description
X063-0005	4.5mm x 30mm Dual Thread Screw
X063-0006	4.5mm x 35mm Dual Thread Screw
X063-0007	4.5mm x 40mm Dual Thread Screw
X063-0008	4.5mm x 45mm Dual Thread Screw
X063-0015	4.5mm x 25mm Single Thread Screw
X063-0016	4.5mm x 30mm Single Thread Screw
X063-0017	4.5mm x 35mm Single Thread Screw
X063-0018	4.5mm x 40mm Single Thread Screw
X063-0019	4.5mm x 45mm Single Thread Screw
X063-0022	10mm Washer

INDICATIONS: See Package Insert for a more complete listing of indications, contraindications, warnings, precautions, and other important information.

LIMITED WARRANTY and DISCLAIMER: Xtant Medical products have a limited warranty against defects and workmanship and materials. Any other express or implied warranties, including warranties of merchantability or fitness, are disclaimed. WARNING: In the USA, this product has labeling limitations. See package insert for complete information. CAUTION: USA Law restricts these devices to sale by or on the order of a physician.

Fixcet is a trademark of Xtant Medical. © Xtant Medical. All Rights Reserved. FM-C-MRK-129 Rev C 11/21



888.886.9354 ✓ cs@xtantmedical.com xtantmedical.com

♀ 664 Cruiser Ln, Belgrade, MT 59714